



Dekoda™ by KOHLER® Health | K-41480

Health tracking device

Features

- Includes Dekoda tracker with battery, magnetic fit pad, authenticator remote and wall-mount bracket, charging block and cable, and quick start guide
- Patented spectroscopy technology and advanced optical sensors observe how light interacts with your waste to analyze gut health, hydration, and the presence of blood
- Pair with the required KOHLER® Health app for daily insights, note dietary habits, track baseline shifts, and share reports directly with your care team
- Fingerprint authentication and secure encrypted data at every step for ultimate privacy
- Use alone or add a family membership for up to six private and separate profiles
- Compatible with most KOHLER White toilets
- Sleek, minimalistic design
- KOHLER Health Membership required (sold separately on Kohlerhealth.com) and auto renewed on a monthly or annual basis until you cancel. Learn more at kohlerhealth.com/terms-of-service



Available Colors/Finishes

Color tiles intended for reference only.

0 White



1-800-4KOHLER (1-800-456-4537)

Kohler Co. reserves the right to make revisions without notice to product specifications.

For the most current specification sheet, go to www.kohler.com USA or www.kohler.ca Canada.

8-06-2026 22:00 - US-CA

THE BOLD LOOK
OF **KOHLER**

Codes/Standards

None Applicable

Technical Information

All product dimensions are nominal.

Height: 3-7/16" (87 mm) - 3-7/16"
(87 mm)

Required Electrical Service

One circuit required.

240 V, 0.5 A, 50/60 Hz

Warranty

See website for detailed warranty information.

KOHLER® One-Year Limited Warranty




1-800-4KOHLER (1-800-456-4537)

Kohler Co. reserves the right to make revisions without notice to product specifications.

For the most current specification sheet, go to www.kohler.com USA or www.kohler.ca Canada.

8-06-2026 22:00 - US-CA



Installation Notes

Install this product according to the installation instructions.

Some toilet models require additional spacers to create the proper gap for Dekoda. You can use seat bumpers or spacers to achieve this spacing by attaching them to the underside of the toilet seat. Please reach out to

support@kohlerhealth.com for support

Please reach out to support@kohlerhealth.com for support



1-800-4KOHLER (1-800-456-4537)

Kohler Co. reserves the right to make revisions without notice to product specifications.

For the most current specification sheet, go to www.kohler.com USA or www.kohler.ca Canada.

8-06-2026 22:00 - US-CA